6 - Precautions for use

The user must have a thorough understanding of devices, instruments and surgical technique before the operation.

Inappropriate selection, positioning or securing of implants can lead to abnormal strain and reduce their useful life. Implantation must always be performed with the specific instruments provided by X.NOV® and in strict adherence with the surgical technique. The cementing technique is very important: comply with the surgical technique provided by X.NOV®. **It is strictly forbidden to make any changes to the implant.**

The procedure must be followed in all cases by suitable post-surgical care and rehabilitation.

Before the operation, check that:

- All the implants that make up the knee prosthesis are available and compatible with each other (check the characteristics indicated on the labelling and leaflet).
- The instruments to be used for the implantation are available and functional. Non-functional or damaged instruments can cause improper positioning of the implants and defectiveness.
- Plan for implants that are at least one size larger and smaller than those determined before the procedure.

General warnings:

- Only remove the implant from the package (in compliance with the usual sterile technique) once the correct size has been determined.
- Pay particular attention to all polished surfaces acting as friction surfaces, and to coated or treated surfaces.
- Never use implants that have been damaged, explanted, subjected to knocks or friction against a hard material.
- Use of cemented versions is recommended in fragile and/or osteoporotic bones.
- Carefully check freedom of movement in order to detect incorrect fixations, instabilities or impingement; correct if necessary. Foreign particles can cause excessive wear-and-tear of the metal/plastic interface. The presence of osteophytes, cement residue and/or bone debris can cause the dislocation of the device or limit movements and make them painful.
- If cement is used, refer to the manufacturer's user's manual.

Attention:
Clean the surfaces well. They must be free from any foreign bodies prior to assembly of the various components and prior to the closure of the surgical site.

DO NOT IMPLANT:
- if the expiry date has been passed
- if the packaging is damaged

7 - Contraindications

- Microcrystalline joint disease.
- Local or systemic, acute or chronic infection.
- Severe muscle, nerve or vessel deficiencies affecting the limb in question.
- Mental or neuromuscular disorders that could constitute an unacceptable risk for the patient and that could cause post-operative complications.
- Insufficient bone stock or poor bone quality, which may affect the stable implantation of the prosthesis
- Any concomitant disease which might influence the implant function.
- Known allergy to implant components.
- Serious diseases that involve a risk of dangerous post-operative complications.
- Osteomyelitis.
- Chronic joint instability (significant instability due to an advanced loss of osteo-cartilaginous structure or a lack of integrity of the lateral ligaments).
- Poor skin coverage or persistent skin lesion on the knee joint.

Infection, sepsis and osteomyelitis are absolute contraindications.

8 - Factors likely to compromise the success of the implantation

- Intense physical activities, risky activities or repetitive movements that expose the implant to excessive stress.
- Allergy to the materials of the implant.
- Medical history of infections and falls.
- BMI > 30 kg/m2.
- Diabetes.
- Serious ligament instability and untreatable muscle contractures.
- Severe osteoporosis.
- Uncooperative or neurologically impaired patient.
- Insufficient bone quality.
- Direct contact of hydroxyapatite-coated implants with acidic local antibiotics.
- Metabolic diseases.
- Drug addiction and/or tendency to abuse drugs, medicinal products, alcohol or tobacco.
- Local bone tumours.
- Flessum > 20°.
- Excessive deformity of the patient's limb (genu varum or severe genu valgum ≥ 30°).
- Damage to implants during surgery.
- Interposition of a foreign body (blood, debris, shavings, etc.) between the implants.
- Pre-operative planning incorrect or not conducted.
- Incorrect selection of the size or type of implants used.
- Non-compliance with the surgical technique.
- Incorrect positioning of the implants.
- Use with an unsuitable instrument or device.
- Major misalignment in the frontal and sagittal plane.

9 - Adverse effects

- Infection.
- Haematoma.
- Haemarthrosis.
- Hyarthrosis.
- Cardiovascular and thrombo-embolic disorders: venous thrombosis, phlebitis and pulmonary or fatty embolism.
- Anaphylactic shock.
- Subchondral tibial collapse.
- Lengthening or shortening of the limb caused by incorrect positioning of the implants.
- Peri-prosthetic fracture and/or fracture of the implant.
- Luxation or subluxation of one or several implants.
- Loosening of one or more components.
- Migration of the implant and loosening of the prosthesis
- Osteolysis.
- Wear-and-tear of components.
- Reaction of tissue to debris from wear-and-tear.
- Sensitivity reactions to the implanted materials.
- Metallosis.
- Soft tissue damage.
- Loss of mobility.
- Pain.
- Neuro-algo-dystrophic syndrome.
- Noises, squeaking.
- Wound, poor healing.
- Infection.
- Revision.

And all of the traditional and known complications of prosthetic surgery.

Repeat surgery may be necessary as a result of the inevitable wear and tear of the device, or any separation that may occur over time.

10 - Storage and handling

Implants must be stored unopened in their original packaging in a cool, dry place away from light. They must not be exposed to ionising radiation, extreme temperatures or particle contamination.

Implants are individually packed in protective and sterile double-barrier packaging. Make sure that the protective packaging is not damaged before removing the implant, as such damage can compromise the sterility of the implant.

Only remove the devices from the package immediately before use.

Disposal of the implant should be carried out in accordance with the method of waste disposal appropriate for the healthcare facility.

Warnings on the package as well as pictograms should be taken into account.

11 - Sterility

The implants are delivered sterile (gamma radiation at a minimum dose of 25 kGy).

They must be kept in their sealed original packaging until use.

Check the sterility expiry date on the label and check the protective packaging for signs of deterioration. If the protective packaging is damaged or the sterility expiry date exceeded, the implant should be sent back to X.NOV®.

The rules for asepsis must be followed when removing the implant from the original packaging and for the entire duration of implantation.

Implants must not be re-sterilised or re-used by the health facility, as this could lead to contamination and/or infection of the patient and damage to the device's properties.

Sterility is compromised if the packaging is damaged.

12 - Patient information

The surgeon must inform the patient of surgical risks, contraindications, possible adverse effects and potential interactions with other procedures or therapeutic or diagnostic devices (e.g.: MRI). He/she must talk to the patient about all the physical and psychological limitations inherent in the use of the implanted device.

Particular attention must be paid to the patient's level of activity, weight, and need for regular medical check-ups.

The surgeon must warn the patient that the implant cannot be as robust, reliable and durable as normal and healthy bones. Therefore, the implant does not allow for the same levels of activity and constraints as those enabled by normal and healthy bones.

The surgeon must explain to the patient that the implant can break or be damaged following a fall or the practice of intense physical activity that could cause trauma or abnormal wear-and-tear to the implant.

If the patient has a profession or activity that uses the implant excessively (e.g. skiing, running, climbing, tennis, body building, etc.), the surgeon must warn the patient and inform him or her that the force exercised on the device and discontinuous and sudden movements that involve the prosthetic joint can cause the failure or reduction in the useful life of the implant or its revision.

Heavy temporary overload in case of falls, accidents or excessive overloads can cause the rupture of the implant, sometime even long after such an event.

The surgeon must tell the patient of the need to keep his/her weight under control and to have regular clinical and radiographic check-ups after the operation.

If the patient experiences pain, loss of mobility or other signs that would lead him or her to believe that the device is not or no longer performing as expected, he or she must consult the practitioner.

13. Compatibility with magnetic resonance imaging (MRI)

Compatibility and safety with MRI environment of X.NOV metallic implants have not been tested to evaluate risk of heating or migration. Risks associated to MRI in presence with a metallic implant include migration and loosening, heating of the implant, heating and damage of the surrounding tissues, artefactual distortion, specifically in the zone surrounding the implant. MRI interpretation may be compromised or made impossible by such artefactual distortions. Patient must be evaluated by a professional knowledgeable of the MRI equipment prior to complete the examination. Patient must inform the professional completing the examination of the presence of his prosthesis in case of MRI examination. Attention to any adverse effect felt by the patient must be maintained all along the MRI examination.

14 - Instrumentation information

The instrumentation for insertion of the implants is supplied in a non-sterile state, unless the contrary is indicated on the packaging. The instrumentation is reusable and must be sterilised by the healthcare facility before use. See the leaflet and the instructions supplied with the instrumentation for information on decontamination and vapour sterilisation of the instrumentation.

Please contact X.NOV® for any additional information on the use of this device.

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Trikompartimentelle posterior stabilisierte Kniegelenkprothese: TRIUM® Femurkomponente / Tibia-Inlay / Tibiasockel / Knie-scheibenimplantat

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