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INSTRUCTIONS FOR USE

Hip prosthesis

CINEOS® femoral head and stem /
5°42' - 12/14 or 6° - 10/12 cone

To be read carefully by the surgeon, theatre staff and all staff involved in the use and handling of the device.

The use of this device is strictly reserved to specialist practitioners and must be carried out in premises that meet the asepsis conditions required and that comply with NF ISO 8828.

Failure to comply with these recommendations for use can reduce the device's useful life, and releases X·NOV® of all liability. Implants must be selected with particular attention to the user, taking account of the relevant indication, the patient's morphology and activity (in particular the type of sport that the patient plays) in order to minimise adverse effects and complications.

X·NOV® bears no liability for complications and/or an abnormally short useful life of the device (these being the sole liability of the user) in the following circumstances:

- incorrect selection of the various implant components,
- modifications made to the implant,
- a poor combination of the various implant components,
- an incorrect indication,
- failure to adhere strictly to the surgical technique recommended by X·NOV®,
- inadequate asepsis,
- failure to comply with the contents of these instructions.

1 - Description

This box contains a femoral stem or a stainless steel or ceramic head CINEOS® intended for the replacement of the natural hip joint in a total hip prosthesis surgery.

This device is an element in the X·NOV® hip prosthesis system and must be used with elements originating from the same manufacturer. This device is intended for single use only.

2 - Material

Metallic femoral head: stainless steel X2CrNiMo18-14-3 in accordance with ISO 5832-1.

Alumina femoral head: NF ISO 6474-1 or NF ISO 6474-2.

Femoral stem: titanium alloy Ti6Al4V in accordance with ISO 5832-3 or X4CrNiMoN21-9-4 in accordance with ISO 5832-9.

Coating of non-cemented stems: hydroxyapatite in accordance with NF ISO 13779-1/2.

3 - Implant selection

Stainless steel femoral heads are available in different embedment lengths on the cone of the femoral stem, which gives rise to different neck lengths. The stainless steel heads are available with an external diameter of 22.2 and 28 mm. For the ceramic heads, the diameters 28, 32, 26 and 40 mm can be used in cups with corresponding diameters. The 5°42' angled cone is impacted on the 12/14 cone of the femoral prosthesis, and the 6° angled cone is impacted on the 10/12 cone.

Femoral stems are available in a graduated range of sizes.

As a rule, the stem size is identical to that of the medullary canal after use of the rasp.

The 5°42'/12-14 cones can take :

- The 5°42'/12-14 stainless steel heads of diameter 22.2 or 28 mm,
- The 5°42'/12-14 alumina heads of diameter 28, 32, 36 or 40 mm.

The 6°/10-12 cones can take:

- The 6°/10-12 stainless steel heads of diameter 22.2 or 28 mm,
- The 6°/10-12 alumina heads of diameter 28, 32, 36 or 40 mm.

The 32, 36 or 40 mm alumina heads must be combined with an XNOV alumina insert with the same diameter.

Implants must be carefully selected during preoperative planning, using the templates provided and the trial implants supplied in the corresponding instrumentation. Poor implant selection can impair the performance of the device and even shorten its useful life.

N.B.: cemented versions are recommended for fragile and/or osteoporotic bones.

According to the current level of knowledge, it is strongly recommended, in the event of revision and/or removal of the prosthesis, that new ceramic components are not implanted in or on prosthetic elements which will remain in place, even if their visual condition seems to be acceptable.

- In case of revision of a broken alumina insert, it is possible to reimplant an alumina insert or a polyethylene insert. If the femoral head must be replaced, the reimplanted head must be made of alumina.
- In case of revision of a broken alumina femoral head, the reimplanted head must be made of alumina.

It is forbidden to re-implant an explanted device.

4 - Indications

- Destruction of the hip joint caused by degenerative or inflammatory conditions.
- Hip fracture
- Avascular necrosis of the femoral head
- Failure or after-effects of previous operations, such as total prosthesis, osteotomy, etc.

5 - Contraindications

- Acute or chronic local or general infections,
- Severe muscle, nerve or vessel deficiencies affecting the limb in question,
- Uncooperative or neurologically impaired patient,
- Major bone defects or poor bone quality, which may affect the stable implantation of the prosthesis,
- Any concomitant disease which might influence the implant function,
- Acidic local antibiotics must not come into direct contact with HA coated implants,
- Allergic reactions,
- Drug addiction and/or tendency to abuse drugs or medication,
- Serious illnesses entailing a risk of dangerous postoperative complications,
- Obesity and overweight.

Infection, septicaemia and osteomyelitis are absolute contraindications.

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6 - Factors likely to compromise the success of the implantation

- Severe osteoporosis,
- Severe deformities, congenital dislocation,
- Local bone tumours,
- Serious ligament instability and untreatable muscle contractures,
- Allergy to implant components,
- Metabolic disorders,
- Medical history of infections and falls,
- Drug dependence,
- BMI > 30kg/m²,
- High-risk activities which could cause implant rupture,
- Physical exertion or repetitive movements that would expose the implant to excessive strain,
- Use with an unsuitable associated device,
- In the event of head revision surgery, use on an existing stem.

7 - Adverse effects

The most common and frequent adverse effects suffered by patients undergoing hip arthroplasty are:

- Migration and loosening of the prosthesis,
 - Dislocation,
 - Infection,
 - Deep vein thrombosis and pulmonary embolism,
 - Infarction,
 - Cardiovascular disorders,
 - Haematoma,
 - Fatigue fracture of the implant,
 - Sensitisation reactions to the materials implanted,
 - Noise, grinding,
 - Wound (poor healing),
 - Component wear and tear,
 - Lengthening or shortening of the limb due to poor positioning,
 - Loosening of one or more components,
 - Rigidity of the hip,
 - Unstable implant,
 - Hydrarthrosis,
 - Hip pain,
 - Osteolysis,
 - Haemarthrosis,
 - Reduced range of movement,
 - Tissue reaction to wear and tear debris.
- Revision surgery may be necessary as a result of the inevitable wear and tear of the device, or any separation that may occur over time.

8 - Handling

Implants must be stored in their original packaging in a cool, dry place away from light. They must not be exposed to ionising radiation, extreme temperatures or particular contamination. Implants are individually packed in double packaging. Packaging must be undamaged. Disposal of the implant should be carried out in accordance with the method of waste disposal appropriate for the healthcare facility.

9 - Precautions for use

Users must have a thorough understanding of devices, instruments and surgical technique before the operation. Inappropriate selection, positioning or securing of implants can lead to abnormal strain and reduce their useful life. Before the operation, check that:

- all the elements making up the hip prosthesis are available,
 - all the elements to be combined are compatible with each other, by checking the properties marked on the labels,
 - all the specific instruments needed to fit the prosthesis are available and in good order.
- Pay particular attention to all polished surfaces acting as friction surfaces, and to coated or treated surfaces. Never use implants that have been damaged, explanted, subjected to knocks or friction against a hard material.

Implantation must always be performed with the specific instruments indicated by the manufacturer.

The protective casing of the head and cone of the prosthesis must not be removed until the final phase of the intervention.

Never wet the ceramic piece. The ceramic heads must not be impacted, but simply positioned with a slight rotating motion.

Good implantation of the femoral head is vital; the success of the arthroplasty and the long-term outcome depend on it.

Before implanting the femoral head, make sure that there is no risk of any particles (organic, cement, etc.) getting between the stem cone and the femoral head.

Before positioning the head on the femoral stem, check the following:

- that the cone of the femoral stem is an X.NOV® cone identical to the cone of the femoral head selected,
 - that the diameter selected is the same as that of the insert,
 - that the stem cone has been carefully cleaned and has no surface defects,
 - that the head is positioned properly in the axis of the male cone.
- Noise/grinding: in literature, up to 20% of cases reported grinding of the hip prostheses with ceramic-on-ceramic interfaces. In order to limit the occurrence of this issue, the implant should be inserted vertically and tilted correctly so that any risk of impingement is avoided. The operation must always be followed by appropriate care and rehabilitation.

DO NOT IMPLANT:
- if the expiry date has been passed,
- if the packaging is damaged.

10 - Sterility

Implants are supplied in a sterile state (ionising radiation minimum 25 kGy).

Implants must not be re-sterilised or re-used, as this could lead to contamination and/or infection of the patient and damage to the device's properties.

Sterility is compromised if the packaging is damaged.

11 - Patient information

The surgeon must inform the patient of the surgical risks, contraindications, possible adverse effects and potential interactions with other therapeutic or diagnostic procedures or equipment, such as MRI. He/she must talk to the patient about all the physical and psychological limitations inherent in the use of the implanted device.

Particular attention must be paid to the patient's level of activity, weight, and need for regular medical check-ups.

The surgeon must warn the patient that the implant cannot be as strong, reliable and long-lasting as normal, healthy bones, and that therefore the patient cannot undertake the same levels of activity and bear the same amount of strain as would be the case with healthy, normal bones. The surgeon must tell the patient that the implant may break or become damaged as a result of intense physical activity or activity that might lead to trauma, and that the implant may need to be replaced in the future.

If the patient has a job or engages in activity that puts excessive strain on the implant (such as skiing, running, climbing, tennis, body-building, etc.), the surgeon must warn the patient and tell him or her that the forces exerted on the device may cause the implant to fail.

The surgeon must tell the patient of the need to keep his/her weight under control and to have regular clinical and radiographic check-ups after the operation.

Patients must consult their doctor if they experience pain (or other symptoms that may indicate that the device is not or no longer performing as it should).

Ceramic implants are sensitive to shock and falls, which are likely to cause ruptures. The literature reports a return rate for hip prostheses with ceramic-on-ceramic interfaces of 1.9% to 3.37% after five years when they are used as intended and in conditions foreseen for this purpose.

12 - Compatibility and safety with MRI environment

Compatibility and safety with MRI environment of X.NOV® metallic implants have not been tested to evaluate risk of heating or migration. Risks associated to MRI in presence of a metallic implant include migration and loosening, heating of the implant, heating and damage of the surrounding tissues, and artefactual distortion, specifically in the area surrounding the implant.

Interpretation of MRI scans may be compromised or made impossible by such artefactual distortion.

Patients must be examined by a professional who is knowledgeable of the MRI equipment prior to undertaking MRI scanning.

Patients must inform the professional responsible for their examination of the presence of their prosthesis in case of MRI scanning. Attention to any adverse effects that may be felt by the patient must be maintained throughout the MRI examination.

13. Instrumentation information

The instrumentation for inserting the implants is supplied in a non-sterile state. It is reusable and must be sterilised by the healthcare facility before use. See the instructions given with the instrumentation kit for information on how to decontaminate and sterilise the instrumentation.

Please contact the manufacturer for all additional information on the use of this device.

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Hüftprothese CINEOS®-Femurkopf und Femurstiel / Konuswinkel 5°42' - 12/14 oder Konuswinkel 6° - 10/12

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